



EFFECTIVE LABORATORY SUPPLIES CO

Balances

GX SERIES

Super Hybrid Sensor, Built-in
Motorized Calibration Weight




- Fast, new Super Hybrid Sensor (SHS) gives you the ultimate in reponse speed – Just ONE Second!
- Resolution Range among models: 1/60,000 to 1/600,000
- Internal Calibration and Internal Calibration Mass makes calibrating your GX easier
- Automatic Self-Calibration calibrates the balance when it detects changes in the ambient temperature
- Large, bright Vacuum Fluorescent Display – see the display even in the dark
- Percentage, Counting Functions & Comparator functions
- ACAI Counting Function for highly accurate counting
- Bright, blue Vacuum Fluorescent Display
- GLP/GMP/ISO Compliance

Model	Weighing Capacity	Maximum Display	Resolution	Repeatability/ Std. Dev.	Linearity
GX-2000	2100g	2100.84 g	0.01 g	0.01 g	0.02g
GX-4000	4100g	4100.84 g	0.01 g	0.01 g	0.02g
GX-6000	6100g	6100.84 g	0.01 g	0.01 g	0.03g
GX-8000	8100g	8100.84 g	0.01 g	0.01 g	0.01g
GX-10000	11000g	11000.84 g	0.01 g	0.01 g	0.01g

Contact

www.effectivelab.co.za • info@effectivelab.co.za • Telephone: (011) 434-1450/1/2 • Fax: (011) 683-1752



EFFECTIVE LABORATORY SUPPLIES CO

Balances

EKi SERIES

High resolution up to 1/60,0000 Resolution (EK-610i, EK-6100i)



- Wide selection of 16 different models 30g ~ 12kg
- Large LCD display (16mm height) with backlight
- Control zero, mode selection and data output with standard RS-232C
- Triple range weighing for EW-i series
- 9 built-in weighing units plus counting & percentage with user priority start-up
- Counting function with Automatic Counting Accuracy Improvement
- Standard comparator function – HI/OK/LO annunciator
- Full Digital Calibration with optional user definable mass setting
- Conformity to GLP with I.D. and serial number outputs
- Auto Power OFF
- Portable-compact size, lightweight, rugged body
- Sealed key panel protection against dust and spills
- Low profile and small footprint for any space requirement
- Overload protection
- Stainless steel weighing pan
- Easy to install rechargeable Ni-MH battery (Optional)
- Optional underhook weighing for EK-4100i/6100i/6000i/12Ki, EW-12Ki
- Leveling feet and level bubble for accurate weighing
- Security ring to prevent theft

Model	Capacity	Resolution	Pan Size
GX-200 ¹	210g	0.001g	128 x 128 mm
GX-400 ¹	410g	0.001g	128 x 128 mm
GX-600 ¹	610g	0.001g	128 x 128 mm
GX-2000	2100g	0.01g	165 x 165 mm
GX-4000	4100g	0.01g	165 x 165 mm
GX-6100	6100g	0.01g	165 x 165 mm
GX-6000	6100g	0.1g	165 x 165 mm
GX-8000	8100g	0.1g	165 x 165 mm

Contact

www.effectivelab.co.za • info@effectivelab.co.za • Telephone: (011) 434-1450/1/2 • Fax: (011) 683-1752



EFFECTIVE LABORATORY SUPPLIES CC

Balances



EJ/EQ/EA

COMPACT BALANCES

- Fast and Stable Weight
- Multi selectable weighing units
Gram, lb, oz, Troy oz, Grain, Carat, Tola, Tael (Taiwan, HK) Momme
- Piece counting and Percentage functions
- Bright and distinct LCD display
- Battery Powered and Auto Off
- Low battery indication
- Leveling feet and level bubble for accurate weighing
- Portable, compact size, light weight, rubbed body
- Stainless steel weighing pan
- Overload protection

Model	Gram	Carat (ct)	Maximum count (pieces)	Minimum unit weight (g)	Sample size (pieces)	Min. Div. (percent)	Min. 100% weight	Pan size	Power
EQ-1200	1200.00x0.01	6000.00x0.05	120,000	0.01	10, 20, 50, 100, 200	0.1%	1g	118mm dia	AC adapter (9V/500mA) or rechargeable battery (6V)
EJ-600	600.00x0.01	3000.00x0.05	60,000	0.01	10, 20, 50, 100, 200	0.1%	1g	118mm dia	
EA-300	300.00x0.01	1500.00x0.05	30,000	0.01	10, 20, 50, 100, 200	0.1%	1g	118mm dia	
EA-600	600.00x0.02	3000.00x0.1	30,000	0.02	10, 20, 50, 100, 200	0.1%	2g	118mm dia	
EA-1500	1500.00x0.05	7500.0x0.2	30,000	0.05	10, 20, 50, 100, 200	0.1%	5g	118mm dia	
EA-3000	3000.0x0.1	15000.0x0.5	30,000	0.1	10, 20, 50, 100, 200	0.1%	10g	118mm dia	

Contact

www.effectivelab.co.za • info@effectivelab.co.za • Telephone: (011) 434-1450/1/2 • Fax: (011) 683-1752