



EFFECTIVE LABORATORY SUPPLIES CC

Filter Presses

ELS Filter Press

- Heavy Duty
- Damax certified stainless steel.
- The tubes for the pots are all certified to a minimum pressure of 20 bar, some four times higher than manufacturer specified pressure utilised.
- The specified max pressure is 5 bar.
- All filter presses come standard with pressure relief valves (max 5 bar).

Filter Press Sizes

100Diam x 300mm High - 2.8Litre
150Diam x 300mm High - 6.2Litre
200Diam x 300mm High - 10.7Litre
250Diam x 300mm High - 16.4Litre
300Diam x 300mm High - 23.4Litre

Available in standard design as well as swing barrel for ease of use – no heavy barrel to lift.

Contact

www.effectivelab.co.za • info@effectivelab.co.za • Telephone: (011) 434-1450/1/2 • Fax: (011) 683-1752

