



EFFECTIVE LABORATORY SUPPLIES CC

Refractometers



Model: PAL-1

- Compact in size allowing for portability
- Measurement is easy. Simply put a drop of a sample on the prism surface. Press START and the Brix value will be displayed in just 3 seconds
- Digital LCD display negates the worry of mis-reading the values
- The design allows the sample to be wiped off easily
- Automatic Temperature Compensation allows measurements to be taken regardless of the ambient temperature. Measurements of samples can be taken during the heating process. (The measured Brix value will stabilize after a few repeated measurements)

International Protection Class: IP 65
Dust-tight and Protected against water jets

Model	Measure. Range	Resolution	Measure. Accuracy	Measure. Temp.	Ambient Temp.	Sample Vol.	Measure. Time	Power	Dimen. (mm)	Weight
PAL-1	Brix 0.0 to 53.0%	Brix 0.1%	Brix \pm 0.2%	10 - 60 °C	10 - 40 °C	0.3ml	3 sec.	2xAAA Batt.	55x31x109	100g

Contact

www.effectivelab.co.za • info@effectivelab.co.za • Telephone: (011) 434-1450/1/2 • Fax: (011) 683-1752