



EFFECTIVE LABORATORY SUPPLIES CC

Trolleys



Stainless Steel Trolleys

This proven range of laboratory and general purpose trolleys are characterized by their robust construction and hygienic ,easy to clean design. Large castor wheels,with vibration absorbing tires,ensure ease of manoeuvrability and no snagging on floor joints or lift entrances. The top trays are at a convenient working height to allow for easy transfer of items to and from any standard bench.

Model No	Number of Trays	Distance between trays	Overall dimensions WxDxH (mm)	Packed Mass (kg)
552	2	510	900 x 560 x 970	20
554	3	200	900 x 560 x 970	25

Features:

- Supplied in knock down form
- Units easy to assemble, instructions provided
- 100mm easy glide ,silent running, heavy duty castors
- Non marking tyres to ensure floors remain clean
- Carrying capacity 50 kg per shelf
- Stainless steel trays
- Sturdy mild steel framework with scratch resistant epoxy powder coated finish

Contact

www.effectivelab.co.za • info@effectivelab.co.za • Telephone: (011) 434-1450/1/2 • Fax: (011) 683-1752



EFFECTIVE LABORATORY SUPPLIES CC

Metal Trolley



- Available in a choice of 2 or 3 shelves
- Standard enamel painted blue
- 100mm castors – 2 fixed and 2 swivel

Model No	Number of Trays	Overall dimensions L x W x H (mm)	Packed Mass (kg)
TSP36/2	2	900 x 535 x 910	250
TSP36/3	3		
TSP42/2	2	1000 x 610 x 1050	250
TSP42/3	3		