



EFFECTIVE LABORATORY SUPPLIES CC

Ultrasonic Cleaners

Model: UMC



- High grade stainless steel insert
- Power supply cable
1.1 metre long cable and plug
- Electromechanical timer
- Temperature controlled heaters


OPTIONAL EXTRAS

- High grade stainless steel baskets
- High grade stainless steel lids
- Special holders to suit customer needs

Model	Nominal Capacity	Inner Tank Dimensions LxWxD (mm)	Overall Dimensions LxWxH (mm)	Effective Liquid Depths	Approx. Mass	Standard Basket Size (mm)	Timer	Heater Controller (Optional)	Power
UMC20	20 liter	495x290x150	515x330x320	110mm	21kg	430x270x80	0 - 30 min	0-40°C	400W
UMC10	10 liter	300x240x160	310x268x303	110mm	14kg	240x210x80	0 - 30 min	0-40°C	200W
UMC5	5 liter	295x150x150	310x175x303	100mm	8kg	130x240x80	0 - 30 min	0-40°C	100W
UMC2	2 liter	150x140x130	180x160x245	105mm	3.3kg	110x110x180	no timer	no heater	50W

Contact

www.effectivelab.co.za • info@effectivelab.co.za • Telephone: (011) 434-1450/1/2 • Fax: (011) 683-1752



EFFECTIVE LABORATORY SUPPLIES CC

Ultrasonic Cleaners

Model: USD300SH-10L


- Stainless steel surround
- Stainless steel lid
- Digitally adjustable timer
- Digitally adjustable heater
- Stainless steel basket included
- 230V supply



Power	Number of Transducers	Tank size	Dimensions LxWxD (mm)	Unit Dimensions LxWxH (mm)	Weight	Timer Range	Heater Range
300W Ultrasonic + 200W heating	6 @ 40kHz	10 liter	295x240x150	327x265x390	8.5kg	0 - 30 min	5 -70°C

Contact

www.effectivelab.co.za • info@effectivelab.co.za • Telephone: (011) 434-1450/1/2 • Fax: (011) 683-1752



EFFECTIVE LABORATORY SUPPLIES CC

Ultrasonic Cleaners

Model: USD500SH-20L



- Stainless steel surround
- Stainless steel lid
- Digitally adjustable timer
- Digitally adjustable heater
- Stainless steel basket included
- 230V supply

Power	Number of Transducers	Tank size	Dimensions LxWxD (mm)	Unit Dimensions LxWxH (mm)	Weight	Timer Range	Heater Range
500W Ultrasonic + 400W heating	10 @ 40kHz	20 liter	495x295x150	530x325x455	16.5kg	0 - 30 min	5 -70°C

Contact

www.effectivelab.co.za • info@effectivelab.co.za • Telephone: (011) 434-1450/1/2 • Fax: (011) 683-1752